

Product Datasheet

Vazyme - VAHTS Universal DNA Library Prep Kit for Illumina V4 (ND610-03)

Description	VAHTS Universal DNA Library Prep Kit for Illumina V4 is a library preparation kit targeted for Illum...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Wide Compatibility: Compatible with FFPE and cfDNA samples of low quality in addition to conventional samples. High Efficiency: Higher library conversion rate. Accurate Detection: Accurate detection of low-frequency mutations.
Expiration Date	12 months from date of receipt.

