

## Product Datasheet

# Vazyme - VAHTS PCR-Free DNA Adapters Set 2 for MGI (NM10902)

<b>Description</b>	VAHTS PCR-Free DNA Adapters for MGI is a specialized kit for MGI high-throughput sequencing platform...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Apply to MGI platform for PCR-Free DNA library preparation. Providing 96 different single indexes.
<b>Expiration Date</b>	12 months from date of receipt.

