

Product Datasheet

Vazyme - VAHTS HiFi Universal Amplification Mix for MGI (SI) (NM618-01)

Description	VAHTS HiFi Universal Amplification Mix for MG1 (SI) is an upgraded version of VAHTS HiFi Amplificati...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Excellent Amplification Uniformity. Extensive Template Compatibility. Excellent Storage Stability.
Expiration Date	12 months from date of receipt.

