

Product Datasheet

Vazyme - VAHTS Library Quantification Kit for Illumina [Out of stock] (NQ101-NQ104)

Description	VAHTS Library Quantification Kit for Illumina is specially designed for accurate quantification of t...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. Master Mix and ROX should be protected from light.
Note	For research use only
Application notes	Contains all components quantified. New dye method qPCR master mix VAHTS SYBR qPCR Master Mix. High amplification efficiency and good GC content compatibility. Accurate quantification of standard products, CV value between batches 5%. Carefully optimized ROX dye concentration, providing ROX premixed version. Compatible with all mainstream quantitative PCR instruments.
Expiration Date	12 months from date of receipt.

