

Product Datasheet

Vazyme - dNTP Mix (10 mM each) (P031-01)

Description	The dNTP Mix is an equimolar mixture of high purity dATP, dTTP, dCTP and dGTP, which can be directly...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	PCR companion: dNTP Mix contains premixed aqueous solutions of dATP, dCTP, dGTP and dTTP, each at a final concentration of 10 mM Quality control: > 99% purity confirmed by HPLC; free from qPCR, PCR, and reverse transcription inhibitors
Expiration Date	12 months from date of receipt.

