

Product Datasheet

Vazyme - dUTP (100 mM) (P033-01)

Description dUTP/dATP/dTTP/dCTP/dGTP are high-purity single component products. The concentration of products wa...

Storage Store at -30 ~ -15°C and transport at ≤0°C

Note For research use only

Application notes PCR companion: Individual components for selection, configurable with different concentrations of dNTPs Quality control: > 99% purity confirmed by HPLC; free from qPCR, PCR, and reverse transcription inhibitors

Expiration Date 12 months from date of receipt.

