

## Product Datasheet

### Vazyme - dGTP (100 mM) (P037-01)

**Description** dUTP/dATP/dTTP/dCTP/dGTP are high-purity single component products. The concentration of products wa...

**Storage** Store at -30 ~ -15°C and transport at ≤0°C

**Note** For research use only

**Application notes** PCR companion: Individual components for selection, configurable with different concentrations of dNTPs Quality control: > 99% purity confirmed by HPLC; free from qPCR, PCR, and reverse transcription inhibitors

**Expiration Date** 12 months from date of receipt.

