

## Product Datasheet

# Vazyme - Taq DNA Polymerase (Mg<sup>2+</sup> plus Buffer) (P101-03)

<b>Description</b>	This product is expressed by Escherichia coli cloned with Thermus aquaticus DNA Polymerase gene, and...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Classic and reliable performance: Optimized buffer system, easy to get ideal amplification results quickly
<b>Expiration Date</b>	12 months from date of receipt.

