

## Product Datasheet

### Vazyme - 2 × Taq Master Mix (P111-03)

<b>Description</b>	This product contains Taq DNA Polymerase, dNTP, and an optimized buffer system. It only needs to add...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Classic and reliable performance: Optimized buffer system, easy to get ideal amplification results quickly Convenient operation: Ready-to-use mix; pre-mixed formulation saves time and reduces contamination
<b>Expiration Date</b>	12 months from date of receipt.

