

Product Datasheet

Vazyme - 2 × Taq Master Mix (Dye Plus) (P112-03)

Description	This product contains Taq DNA Polymerase, dNTP, and an optimized buffer system. It only needs to add...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Classic and reliable performance: Optimized buffer system, easy to get ideal amplification results quickly Convenient operation: Include tracking dye, direct loading of PCR products on a gel
Expiration Date	12 months from date of receipt.

