

Product Datasheet

Vazyme - Champagne Taq antibody [Out of stock] (P121-01)

Description	Champagne Taq antibody is a monoclonal antibody against Taq DNA polymerase, which can inhibit the po...
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Reliable performance: Thermostable Taq monoclonal antibody; high specificity Quality Control: Suspension cell fermentation production; stringent quality control Good compatibility: Compatible with different brands of Taq enzymes
Expiration Date	12 months from date of receipt.

