

Product Datasheet

Vazyme - 2 × Vazyme LAMP Master Mix (Dye Plus) (P312-03)

Description	Vazyme LAMP DNA Polymerase is a mixed enzyme of Taq DNA Polymerase and a protein containing 3'→5' ...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Good performance: ~6X higher than Taq DNA Polymerase; higher amplification efficiency Convenient operation: Ready-to-use mix; include tracking dye, direct loading of PCR products on a gel Additional features: 3'→5' exonuclease activity Broad Compatibility: Amplify fragments up to 21 kb (gDNA)
Expiration Date	12 months from date of receipt.

