

## Product Datasheet

# Vazyme - AceTaq DNA Polymerase [Out of stock] (P401-MD1)

**Description** AceTaq DNA Polymerase is a chemically modified Taq DNA Polymerase that is completely blocked at room...

**Storage** Store at -30 ~ -15°C and transport at ≤0°C.

**Note** For research use only

**Application notes** High stability: Chemical modification; complete enzymatic activity below 80°C  
Rapid activation: Activity can be released by heating at 95°C for 5 minutes  
High sensitivity: Suitable for amplifying low-copy genes from complex templates

**Expiration Date** 12 months from date of receipt.

