

## Product Datasheet

# Vazyme - Phanta Super-Fidelity DNA Polymerase (P501-d1)

<b>Description</b>	Phanta Super-Fidelity DNA Polymerase is a new generation of super-fidelity DNA polymerase modified b...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C. ▲ Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	High fidelity: ~52X higher than Taq DNA Polymerase Fast amplification: Extension speed up to 15 - 30 sec/kb Broad Compatibility: Suitable for various types of templates Convenient operation: Robust amplification with minimal optimization
<b>Expiration Date</b>	12 months from date of receipt.

