

Product Datasheet

Vazyme - Phanta Super-Fidelity DNA Polymerase (P501-d2)

Description	Phanta Super-Fidelity DNA Polymerase is a new generation of super-fidelity DNA polymerase modified b...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. ▲ Avoid repeated freezing and thawing.
Note	For research use only
Application notes	High fidelity: ~52X higher than Taq DNA Polymerase Fast amplification: Extension speed up to 15 - 30 sec/kb Broad Compatibility: Suitable for various types of templates Convenient operation: Robust amplification with minimal optimization
Expiration Date	12 months from date of receipt.

