

## Product Datasheet

# Vazyme - Phanta Super-Fidelity DNA Polymerase (P501-d3)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Phanta Super-Fidelity DNA Polymerase is a new generation of super-fidelity DNA polymerase modified b...  |
| <b>Storage</b>           | Store at -30 ~ -15°C and transport at ≤0°C. ▲ Avoid repeated freezing and thawing.   |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | High fidelity: ~52X higher than Taq DNA Polymerase<br>Fast amplification: Extension speed up to 15 - 30 sec/kb<br>Broad Compatibility: Suitable for various types of templates<br>Convenient operation: Robust amplification with minimal optimization |
| <b>Expiration Date</b>   | 12 months from date of receipt.  |

