

Product Datasheet

Vazyme - Phanta Max Super-Fidelity DNA Polymerase (P505-d1)

Description	Phanta Max Super-Fidelity DNA Polymerase is an upgraded version of Phanta Super-Fidelity DNA Polymerase.
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Ultra-high fidelity: ~128X higher than Taq DNA Polymerase Longer amplicon lengths: Up to 40 kb simple templates, 10 kb cDNA and 20 kb gDNA Robust performance: Superior specificity, sensitivity, and yields Broad Compatibility: Suitable for complex samples & GC-rich systems
Expiration Date	12 months from date of receipt.

