

## Product Datasheet

### Vazyme - 2 × Phanta Flash Master Mix (P510-02)

<b>Description</b>	2 × Phanta Flash Master Mix is a new generation superior enzyme based on Phanta Flash Super-Fidelit...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	High fidelity: ~81X higher than Taq DNA Polymerase Flash amplification: Extension speed up to 5 sec/kb Robust performance: Superior specificity, sensitivity, and yields Broad compatibility: Suitable for complex samples & GC-rich systems
<b>Expiration Date</b>	12 months from date of receipt.

