

Product Datasheet

Vazyme - 2 × Phanta Flash Master Mix (Dye Plus (P520-03))

Description	2 × Phanta Flash Master Mix is a new generation superior enzyme based on Phanta Flash Super-Fidelit...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	High fidelity: ~81X higher than Taq DNA Polymerase Flash amplification: Extension speed up to 5 sec/kb Robust performance: Superior specificity, sensitivity, and yields Broad Compatibility: Suitable for complex samples & GC-rich systems
Expiration Date	12 months from date of receipt.

