

## Product Datasheet

# Vazyme - 2 × Phanta Max Master Mix (Dye Plus (P525-01))

<b>Description</b>	Phanta Max has added unique elongation factors, specificity promoting factors and plateau uninhibiti...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Ultra-high fidelity: ~128X higher than Taq DNA Polymerase Longer amplicon lengths: Up to 40 kb simple templates, 10 kb cDNA and 20 kb gDNA Robust performance: Superior specificity, sensitivity, and yields Broad Compatibility: Suitable for complex samples & GC-rich systems
<b>Expiration Date</b>	12 months from date of receipt.

