

Product Datasheet

Vazyme - Single Cell Sequence Specific Amplification Kit [Out of stock] (P621-01)

Description	The Single Cell Sequence Specific Amplification Kit is a one-step RT-PCR-based amplification method....
Storage	Store at -20°C for up to 18 months.
Note	For research use only
Application notes	Suitable for 1-10,000 cells; One-step RT-PCR pre-amplification; Up to 500 target genes can be amplified; The amplified product is suitable for analysis by any real-time fluorescent quantitative PCR instrument; High sensitivity, low cost
Expiration Date	Please enquire.

