

Product Datasheet

Radox - Toltrazuril (PAS9456)



Description	Anti-Toltrazuril Antibody
Species/Host	Sheep
Conjugation	Unconjugated
Immunogen	Toltrazuril-BTG
Target	Toltrazuril
Preservatives	0.09% Sodium Azide.. 20mM Phosphate, 150mM Sodium Chloride. pH 7.291
Form/Appearance	Clear Liquid
Concentration	Lot dependent (Determined @ 280nm.)
Storage	Can be stored for up to 3 months at +2°C - +8°C. For long term storage, aliquot and store at ≤-20°C. Avoid repeated freeze/thaw cycles. Product should be protected from light exposure.
Note	For research use only
Application notes	The antibody is suitable for the development of immunoassays or immunoaffinity purification columns
Isotype	IgG
Clonality	Polyclonal
Sensitivity	10 ng/mL Toltrazuril produces 32% inhibition in a competitive ELISA employing Toltrazuril polyclonal antibody.
Dilution Range	0.625 µg/mL
Expiration Date	12 months from date of receipt.