

Product Datasheet

Radox - Zeranol (PAS9464)

Description	Anti-Zeranol Antibody
Species/Host	Sheep
Conjugation	Unconjugated
Immunogen	Zearalanone-BTG
Target	Zeranol
Preservatives	0.09% Sodium Azide.. 20mM Phosphate, 150mM Sodium Chloride. pH 7.317
Form/Appearance	Clear Liquid
Concentration	Lot dependent (Determined @ 280nm.)
Storage	Can be stored for up to 3 months at +2°C - +8°C. For long term storage, aliquot and store at ≤-20°C. Avoid repeated freeze/thaw cycles. Product should be protected from light exposure.
Note	For research use only
Application notes	The antibody is suitable for the development of immunoassays or immunoaffinity purification columns
Isotype	IgG
Clonality	Polyclonal
Sensitivity	10 ng/mL Zeranol produces 82% inhibition in a competitive ELISA employing Zeranol polyclonal antibody.
Dilution Range	2.5 µg/mL
Expiration Date	12 months from date of receipt.

