

Product Datasheet

Radox - Medroxyprogesterone Acetate (PAS9538)

Description	Anti-Medroxyprogesterone Acetate Antibody
Species/Host	Sheep
Conjugation	Unconjugated
Immunogen	Medroxyprogesterone Acetate-BTG
Target	Medroxyprogesterone Acetate
Preservatives	0.09% Sodium Azide.. 20mM Phosphate, 150mM Sodium Chloride. pH 7.57
Form/Appearance	Clear Liquid
Concentration	Lot dependent (Determined @ 280nm.)
Storage	Can be stored for up to 3 months at +2°C - +8°C. For long term storage, aliquot and store at ≤-20°C. Avoid repeated freeze/thaw cycles. Product should be protected from light exposure.
Note	For research use only
Application notes	The antibody is suitable for the development of immunoassays or immunoaffinity purification columns
Isotype	IgG
Clonality	Polyclonal
Sensitivity	10 ng/mL Medroxyprogesterone Acetate produces 99% inhibition in a competitive ELISA employing Medroxyprogesterone Acetate polyclonal antibody.
Dilution Range	0.0938 µg/mL
Expiration Date	12 months from date of receipt.

