

## Product Datasheet

### Vazyme - AceQ qPCR Probe Master Mix (Q112-03)

<b>Description</b>	The AceQ qPCR Probe Master Mix is specially designed for probe-based qPCR. The core component Ace Ta...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C
<b>Note</b>	For research use only
<b>Application notes</b>	<a href="https://www.vazymebiotech.com/uploads/Q112-Components.png">https://www.vazymebiotech.com/uploads/Q112-Components.png</a>
<b>Expiration Date</b>	12 months from date of receipt.

