

Product Datasheet

Vazyme - HiScript II One Step qRT-PCR Probe Kit [Out of stock] (Q222-01)

Description HiScript II One Step qRT-PCR Probe Kit is specially designed for qPCR detection using RNA (e.g., RNA...

Storage Store at -30 ~ -15°C and transport at ≤0°C.

Note For research use only

Application notes Ultra-high detection sensitivity: single-digit copies of templates or 0.1 pg total RNA can be detected; Multiple detection: support multiple probe detection

Expiration Date Please enquire.

