

## Product Datasheet

# Vazyme - HiScript III One Step qRT-PCR Probe Kit (Q225-EN01)

**Description** HiScript III One Step qRT-PCR Probe Kit is specially designed for qPCR that directly use RNA (e.g. v...

**Storage** Store at -30 ~ -15°C and transport at ≤ 0°C.

**Note** For research use only

**Application notes** Higher detection sensitivity  
Balanced amplification at high and low template concentrations  
Applicable to multiplex qPCR  
Better storage stability  
Compatible with fast mode

**Expiration Date** 12 months from date of receipt.

