

Product Datasheet

Vazyme - HiScript III One Step qRT-PCR Probe Kit (Q225-EN02)

Description HiScript III One Step qRT-PCR Probe Kit is specially designed for qPCR that directly use RNA (e.g. v...

Storage Store at -30 ~ -15°C and transport at ≤ 0°C.

Note For research use only

Application notes Higher detection sensitivity
Balanced amplification at high and low template concentrations
Applicable to multiplex qPCR
Better storage stability
Compatible with fast mode

Expiration Date 12 months from date of receipt.

