

## Product Datasheet

# Vazyme - HiScript III One Step qRT-PCR Probe Kit (Q225-EN03)

<b>Description</b>	HiScript III One Step qRT-PCR Probe Kit is specially designed for qPCR that directly use RNA (e.g. v...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤ 0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Higher detection sensitivity Balanced amplification at high and low template concentrations Applicable to multiplex qPCR Better storage stability Compatible with fast mode
<b>Expiration Date</b>	12 months from date of receipt.

