

Product Datasheet

Vazyme - AccurSTART One Step RT-qPCR Probe Kit (FOR FAST) (Q231-EN01)

Description	AccurSTART One Step RT-qPCR Probe Kit (FOR FAST) is a probe-based RT-qPCR fast detection kit suitabl...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Rapid amplification Sensitive and efficient Compatible and impurity-tolerable Co-detection of DNA and RNA
Expiration Date	12 months from date of receipt.

