

## Product Datasheet

# Vazyme - AceQ Universal SYBR qPCR Master Mix (Q511-02)

<b>Description</b>	AceQ Universal SYBR qPCR Master Mix is a special premix for qPCR detection using SYBR Green I chimer...
<b>Storage</b>	Store at -30 ~ -15°C and protect from light. Transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Universal ROX Passive Reference Dye: It is suitable for all qPCR instruments. The concentration of ROX does not need to be adjusted on different instruments. Rigorous hot-start DNA polymerase: Novel chemically modified hot-start Ace Taq DNA polymerase with strong amplification performance and high specificity. Optimized reaction system: Inhibit non-specific amplification and primer dimerization to a large extent. Stable well-to-well repeatability: Meet the requirements of high-class magazines.
<b>Expiration Date</b>	12 months from date of receipt.

