

Product Datasheet

Vazyme - AceQ Universal Probe Master Mix V2 (Q513-EN02)

Description AceQ Universal Probe Master Mix V2 is a probe-based reagent for qPCR. The core component AceTaq DNA ...



Storage Store at -30 ~ -15°C and transport at ≤0°C. Keep away from light

Note For research use only

Application notes Superior amplification sensitivity: Hot-start Ace Taq DNA Polymerase provides superior amplification performance and sensitivity. Excellent linear relationship: The amplification efficiency can reach to 99%. Reliable quantification results: It contains a unique anti-contamination system to elimination deviation caused by aerosols. Universal ROX Reference Dye: It contains a special ROX reference dye, which is applicable for almost all qPCR instruments.

Expiration Date 12 months from date of receipt.