

Product Datasheet

Vazyme - HiScript III One Step qRT-PCR Probe 5 × Master Mix (Q611-EN02)

Description	HiScript III One Step qRT-PCR Probe 5 × Master Mix is single-tube RT-qPCR premix that is suitable f...
Storage	Store at -30 ~ -15°C and transport at ≤ 0°C.
Note	For research use only
Application notes	Excellent stability, easy to add samples (one tube mix) Excellent balance of high and low concentration template amplification Support full premix (except probes and primers)
Expiration Date	12 months from date of receipt.

