

## Product Datasheet

# Vazyme - HiScript III One Step qRT-PCR Probe 5 × Master Mix (Q611-EN03)

<b>Description</b>	HiScript III One Step qRT-PCR Probe 5 × Master Mix is single-tube RT-qPCR premix that is suitable f...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤ 0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Excellent stability, easy to add samples (one tube mix) Excellent balance of high and low concentration template amplification Support full premix (except probes and primers)
<b>Expiration Date</b>	12 months from date of receipt.

