

Product Datasheet

Vazyme - AccurSTART One Step RT-qPCR Super PreMix (ONE TUBE) (Q621-EN02)

Description	AccurSTART One Step RT-qPCR Super PreMix (ONE TUBE) is a one-tube probe-based RT-qPCR master mix sui...
Storage	Store at -30 ~ -15°C and transport at ≤ 0°C.
Note	For research use only
Application notes	Primer and probe premixed: Probe and primers can be added to the premix in advance without changing the amplification performance DNA/RNA co-detection: Stable and sensitive amplification of DNA and RNA templates Multiplex amplification: Compatible with 4 - 5 multiplex amplification Sensitive and efficient: Achieve stable detection of single-digit copy templates
Expiration Date	12 months from date of receipt.

