

Product Datasheet

Vazyme - ChamQ Geno-SNP Probe Master Mix (Q811-03)

Description	ChamQ Geno-SNP Probe Master Mix is specially designed for single nucleotide polymorphism (SNP) typin...
Storage	Store at -30 ~ -15°C and transport at ≤0°C. Keep away from light.
Note	For research use only
Application notes	Compatible with direct amplification of blood lysate: SNP typing can be performed using only blood lysate as th template. Broad template input compatibility: It is compatible with template input ranging from 1 - 10 ng. Great GC-rich template compatibility: It accurately genotypes at SNP sites with a GC content of 25 - 73%.
Expiration Date	12 months from date of receipt.

