

## Product Datasheet

# Vazyme - AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) (QL227-EN01)

<b>Description</b>	AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) is a glycerol-free reagent using RNA as a temp...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤ 0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Sensitive and efficient: Excellent amplification performance Excellent tolerance to impurities: Compatible with different samples (e.g. feces, tissue, and swab) Strong lyophilizable stability: Superior stability of accelerated test and freeze-thaw cycles
<b>Expiration Date</b>	12 months from date of receipt.

