

## Product Datasheet

# Vazyme - AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) (QL227-EN03)

**Description** AccurSTART One Step RT-qPCR Probe Kit (Glycerol-free) is a glycerol-free reagent using RNA as a temp...

**Storage** Store at -30 ~ -15°C and transport at ≤ 0°C.

**Note** For research use only

**Application notes** Sensitive and efficient: Excellent amplification performance Excellent tolerance to impurities: Compatible with different samples (e.g. feces, tissue, and swab) Strong lyophilizable stability: Superior stability of accelerated test and freeze-thaw cycles

**Expiration Date** 12 months from date of receipt.

