

Product Datasheet

Vazyme - Taq Pro Multiple Probe qPCR Mix (QN213-EN02)

Description	Taq Pro Multiple Probe qPCR Mix is a master mix for probe qPCR to detect DNA templates (such as DNA ...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Great low-copy templates amplification sensitivity and 'S-shaped' amplification curve Compatible with Fast Program Compatible with multiplex detection Strong anti-contamination ability
Expiration Date	12 months from date of receipt.

