

## Product Datasheet

# Vazyme - Animal Detection Probe Master Mix (QV110-EN01)

<b>Description</b>	Animal Detection Probe Master Mix is a specialized master mix for probe qPCR to detect African swine...
<b>Storage</b>	Store at -30 ~ -15°C and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Superior sensitivity: Novel antibody modified hot-start Taq DNA polymerase with the optimized buffer, increasing the performance of low-concentration template detection. Anti-contamination system: Efficient degradation of aerosol contamination, ensuring reliable qPCR results. Support Fast Program: The detection can be finished within 40 min. Strong storage stability: The performance remains stable after acceleration at 37°C for 9 days and 30 freeze-thaw cycles.
<b>Expiration Date</b>	12 months from date of receipt.

