

Product Datasheet

Vazyme - Animal Detection Probe Master Mix (QV110-EN02)

Description	Animal Detection Probe Master Mix is a specialized master mix for probe qPCR to detect African swine...
Storage	Store at -30 ~ -15°C and transport at ≤0°C.
Note	For research use only
Application notes	Superior sensitivity: Novel antibody modified hot-start Taq DNA polymerase with the optimized buffer, increasing the performance of low-concentration template detection. Anti-contamination system: Efficient degradation of aerosol contamination, ensuring reliable qPCR results. Support Fast Program: The detection can be finished within 40 min. Strong storage stability: The performance remains stable after acceleration at 37°C for 9 days and 30 freeze-thaw cycles.
Expiration Date	12 months from date of receipt.

