

## Product Datasheet

# Vazyme - HiScript II Q RT SuperMix for qPCR [Out of stock] (R222-01)

<b>Description</b>	HiScript II Reverse Transcriptase is a novel reverse transcriptase obtained by in vitro molecular ev...
<b>Storage</b>	Store at -30 ~ -15°C for up to 12 months, and transport at ≤0°C.
<b>Note</b>	For research use only
<b>Application notes</b>	Based on the efficient HiScript II Reverse Transcriptase, the reverse transcription efficiency is higher; Anchored Oligo (dT)23 VN is uniquely designed to ensure reaction specificity and improve the efficiency and success rate of first-strand cDNA synthesis; The ready-to-use SuperMix is easy to use and time-saving
<b>Expiration Date</b>	12 months from date of receipt.

