

Product Datasheet

Vazyme - HiScript-TS 2 × PCR Mix [Out of stock] (RA103-01)

Description	HiScript-TS 2 × PCR Mix is a dedicated reagent for HiScript-TS 5'/3' RACE Kit (Vazyme #RA101), suit...
Storage	Store at -30 ~ -15°C and transport at ≤0°C
Note	For research use only
Application notes	HiScript-TS 5'/3' RACE Kit (Vazyme #RA101) dedicated cDNA ends rapid amplification reagent
Expiration Date	12 months from date of receipt.

