

---

## Product Datasheet

**Radox - Anti Prolactin sdAb clone 9.24E10.F3.H12  
(RAF9651)**



<b>Description</b>	Recombinant single domain antibody expressed in an E. Col...
<b>Species/Host</b>	Camelus
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Human
<b>Target</b>	Prolactin
<b>Preservatives</b>	1x PBS containing 0.09% sodium azide preservative.
<b>Concentration</b>	Typically > 1mg/ml
<b>Storage</b>	These fragments are stable at 4°C. It is recommended that for storage over extended periods they are kept at -20°C and should not be subject to repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	These fragments contain His and c-Myc fusion tags which may be used for detection or immobilisation. <sup>2</sup> Recombinant antibody fragments are suitable for use in ELISA immunoassays, biosensor applications, western blots, immunohistochemistry, flow cytometry, immunoaffinity purification and most other immunological methods*
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 90% assessed by SDS-PAGE.
<b>MW</b>	Approximately 35 kDa
<b>Dilution Range</b>	To be determined by end-user.
<b>Expiration Date</b>	12 months from date of receipt.