

---

## Product Datasheet

**Radox - Anti PD1 (Programmed cell death protein 1) sdAb clone 6F7 (RAF9747)**



**Description**

Recombinant single domain antibody

and expressed in an E. Col...

**Species/Host**

Camelus

**Conjugation**

Unconjugated

**Immunogen**

Human

**Target**

Programmed Cell Death Protein 1 (PD1)

**Preservatives**

1x PBS containing 0.09% sodium azide preservative.

**Concentration**

Typically > 1mg/ml

**Storage**

These fragments are stable at 4°C. It is recommended that for storage over extended periods they are kept at -20°C and should not be subject to repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

These fragments contain His and c-Myc fusion tags which may be used for detection or immobilisation.<sup>2</sup> Recombinant antibody fragments are suitable for use in ELISA immunoassays, biosensor applications, western blots, immunohistochemistry, flow cytometry, immunoaffinity purification and most other immunological methods\*

**Clonality**

Recombinant

**Purity**

> 90% assessed by SDS-PAGE.

**MW**

Approximately 35 kDa

**Dilution Range**

To be determined by end-user.

**Expiration Date**

12 months from date of receipt.