
Product Datasheet

Radox - Anti CTLA4 (Cytotoxic T-Lymphocyte-associated Antigen 4) sdAb clone 8C7 (RAF9763)



Description

Recombinant single domain antibody

Alpaca and expressed in an E. Col...

Species/Host

Camelus

Conjugation

Unconjugated

Immunogen

Human

Target

Cytotoxic T-Lymphocyte-associated Antigen 4 (CTLA4)

Preservatives

1x PBS containing 0.09% sodium azide preservative.

Concentration

Typically > 1mg/ml

Storage

These fragments are stable at 4°C. It is recommended that for storage over extended periods they are kept at -20°C and should not be subject to repeated freeze-thaw cycles.

Note

For research use only

Application notes

These fragments contain His and c-Myc fusion tags which may be used for detection or immobilisation.² Recombinant antibody fragments are suitable for use in ELISA immunoassays, biosensor applications, western blots, immunohistochemistry, flow cytometry, immunoaffinity purification and most other immunological methods*

Clonality

Recombinant

Purity

> 90% assessed by SDS-PAGE.

MW

Approximately 35 kDa

Dilution Range

To be determined by end-user.

Expiration Date

12 months from date of receipt.