
Product Datasheet

Radox - Anti CTLA4 (Cytotoxic T-Lymphocyte-associated Antigen 4) sdAb clone 8G1 (RAF9771)



Description

Recombinant single domain antibody

Species/Host	Alpaca and expressed in an E. Col...
Conjugation	Unconjugated
Immunogen	Human
Target	Cytotoxic T-Lymphocyte-associated Antigen 4 (CTLA4)
Preservatives	1x PBS containing 0.09% sodium azide preservative.
Concentration	Typically > 1mg/ml
Storage	These fragments are stable at 4°C. It is recommended that for storage over extended periods they are kept at -20°C and should not be subject to repeated freeze-thaw cycles.
Note	For research use only
Application notes	These fragments contain His and c-Myc fusion tags which may be used for detection or immobilisation. ² Recombinant antibody fragments are suitable for use in ELISA immunoassays, biosensor applications, western blots, immunohistochemistry, flow cytometry, immunoaffinity purification and most other immunological methods*
Clonality	Recombinant
Purity	> 90% assessed by SDS-PAGE.
MW	Approximately 35 kDa
Dilution Range	To be determined by end-user.
Expiration Date	12 months from date of receipt.