

Product Datasheet

Vazyme - FastPure Viral DNA/RNA Mini Kit Pro (RC323-01)

Description	FastPure Viral DNA/RNA Mini Kit Pro is a column extraction kit, which can extract viral DNA/RNA from...
Storage	Box 1: Store at -30 ~ -15°C and transport at ≤ 0°C. Box 2: Store at 15 ~ 25°C and transport at room temperature.
Note	For research use only
Application notes	Extreme performance: Easily capture trace nucleic acids of virus Wide sample compatibility: Suitable for whole blood, serum, plasma, swab, tissue, alveolar lavage fluid, cell culture supernatant and other samples Efficient and fast: Operate at room temperature for only 12 min
Expiration Date	12 months from date of receipt.

