

## Product Datasheet

# Vazyme - FastPure Universal Plant Total RNA Isolation Kit (RC411-01)

<b>Description</b>	FastPure Universal Plant Total RNA Isolation Kit can be used to quickly extract total RNA from plant...
<b>Storage</b>	Store at 15 ~ 25°C and transport at room temperature.
<b>Note</b>	For research use only
<b>Application notes</b>	Simple 11 min procedure No need for phenol or chloroform and $\beta$ -mercaptoethanol Compatible with plants rich in polysaccharide and polyphenol
<b>Expiration Date</b>	Please enquire.

