

Product Datasheet

Radox - Interleukin-12p40 (IL-12p40) (RCP9122)

Description	Active IL-12 is a p70 disulphide-linked dimer composed of p35 and p40 subunits. The protein is a ple...
Preservatives	25 mM Sodium phosphate pH 7, 50% Glycerol
Form/Appearance	Liquid
Storage	at -20 °C
Tag	N terminal containing a HIS tag
Note	For research use only
Application notes	ELISA, Western Blotting, etc.
Purity	> 95 %, by SDS-PAGE under reducing conditions
Source	E.coli -derived human recombinant protein
MW	39.2 kDa
Uniprot ID	P29460
Expiration Date	6 months from date of receipt.

